



- Double walled, Exterior made of CRCA steel - powder coat finished.
- Inner Leak proof tank made of Stainless Steel – 304 grade (Argon welded).
- High density glass wool insulation to prevent heat loss & ensure constant temperature.
- Pump cum stirring assembly for high accuracy and temperature uniformity.
- Temperature controlled by Dual display Microprocessor based PID Temp. Controller with PT-100 temperature sensor.
- Temperature Range: 5°C above ambient to 99° C.
- Temperature Accuracy: $\pm 0.5^\circ$ C.
- Drainage at the side of the bath is provided as outlet.
- Operates on 230 Volts AC Single phase 50 Hz.

MODEL	WORKING SIZE		RATING
	INCHES	CM	
LCTB-1	12 X 10 X 8	30 X 25 X 20	1500 W.
LCTB-2	18 X 12 X 10	45 X 30 X 25	2500 W.
LCTB-3	18 X 14 X 12	45 X 35 X 30	3000 W.

OPTIONAL ACCESSORIES :

- A. Digital Timer with Temp. Controller of 999 minutes.
- B. External Recirculation.
- C. Glass window.