

HPLC COLUMN STORAGE CABINETS



CSC 30SS

Specification

■ CSC 30

- Capacity : 30 HPLC column
- Dimensions (WxHxD) : 295 x 380 x 400 mm
- No. of Drawers : 5 drawers

■ CSC 60

- Capacity : 60 HPLC column
- Dimensions (WxHxD) : 295 x 760 x 400 mm
- No. of Drawers : 10 drawers

■ CSC 90

- Capacity : 100 HPLC column
- Dimensions (WxHxD) : 295 x 1140 x 400 mm
- No. of Drawers : 15 drawers

■ CSC 120

- Capacity : 120 HPLC column
- Dimensions (WxHxD) : 590 x 760 x 400 mm
- No. of Drawers : 20 drawers

Available in 4 Sizes (Standard)

- Column storage 30, 60, 90 & 120
- Also available 250, 500 & 1000 columns

Stackable

Unique stacking kit allows increase the number of columns you can store at any time

100% Steel Cabinet

Quality, Solid Steel Construction with chrome plated O-Ring handles

Versatile

Store upto 12 x 12.5 cm columns in a single drawer, or a combination of short columns, guard columns & analytical columns. Holds all columns upto 30 cm long

Special size available as per requirement

Stainless Steel (SS) cabinet also available

Colour Options (Red / Blue / Gray / Off White)



CSC 30



CSC 60



CSC 90



CSC 120



CSC 100 / CSC 200 / CSC 250
CSC 500 / CSC 1000