



- Oil baths are ideal for reactions, extraction, analysis & other lab tests.
- Double walled in construction.
- Leak proof argon welded Inner made of Stainless Steel – 304 grade & Exterior made of CRCA Steel -powder coat finish.
- Inner tank of Stainless Steel – 304 / 316 grade & Body made of Stainless Steel – 304 grade Matt finished.
- Walls contain High density cerawool insulation to minimize the heat loss.
- Temperature controlled by Capillary Thermostat with +/- 2° C accuracy.
- Temp. Range 5° C Above ambient to 250° C.
- Supplied without lid.
- Operates on 230 Volts AC Single phase 50 Hz.

**OPTIONAL ACCESSORIES: -**

A)Dual display Microprocessor based PID Temperature Controller with Temp. Accuracy : ± 2 C OR Better.

MODEL	SIZE IN INCH Ø X Ht	CAPACITY	RATING
LHOB-1 S/G	7" X 5"	1 LITERS FLASK	1000
LHOB-2 S/G	8" X 6"	2 LITERS FLASK	1000
LHOB-3 S/G	9" X 7"	3 LITERS FLASK	1500
LHOB-4 S/G	11" X 9"	5 LITERS FLASK	2000
LHOB-5 S/G	13" X 11"	10 LITERS FLASK	2000
LHOB-6 S/G	15" X 13"	20 LITERS FLASK	2500