



- Electric Bunsen Burners are designed & used for evaporation of solutions in china dish & can be used for any routine heating application in the laboratory.
- Heating element is fixed within the conically shaped refractory cup.
- Cup is enclosed in a stainless steel casing and well insulated.
- S. S. casing is mounted on a stainless steel body with perforations to prevent temperature rise.
- Temp. controlled by Energy regulator.
- Operates on 230 Volts AC Single phase 50 Hz.

MODEL	RATING	MAXIMUM TEMP.
LBNB-1	300 W.	700° C

**OPTIONAL ACCESSORIES :**

A. Dual display Microprocessor based PID temperature controller & PT – 100 sensor.