

**OPTIONAL ACCESSORIES :**

A. Dual display Microprocessor based PID temperature controller & PT – 100 sensor.

- Non rusting Aluminum housing - powder coat finish.
- Thick hand knitted fiber glass netting.
- High density glass wool insulation.
- Energy regulator – power controller.
- Maximum surface temperature 350° C.
- Booster heaters for bigger sizes beyond 10 Ltrs.
- High grade chrome plated Nichrome heating elements.
- Operates on 230 Volts AC Single phase 50 Hz.

MODEL	FLASK CAPACITY	WATTS
LHTM-1	100 ml	60
LHTM-2	250 ml	150
LHTM-3	500 ml	200
LHTM-4	1 Litre	300
LHTM-5	2 Litres	450
LHTM-6	3 Litres	500
LHTM-7	5 Litres	600
LHTM-8	10 Litres	2 x 600
LHTM-9	20 Litres	3 x 600