



- Plates are in Rectangular/Round Shape.
- Top plate made form Cast Iron with body made of thick CRCA Steel, powder coat finished.
- Body & Top plate made of Stainless Steel – 304 grade in GMP model.
- Long Lasting Nichrome Elements embedded with porcelain insulation beads.
- Temperature controlled by Thermostat & Heater Switch OR Energy Regulator
- Max. surface temp. 300° C.
- Operates on 230 Volts AC Single phase 50 Hz .

MODEL		PLATE SIZE	WATTS.
ROUND	LHTP-1 S/G	8" Ø	1200
	LHTP-2 S/G	9" Ø	1500
	LHTP-3 S/G	12" Ø	2000

MODEL		PLATE SIZE	WATTS.
RECTANGULAR	LHTP-4 S/G	10" X 12"	1200
	LHTP-5 S/G	10" X 16"	1500
	LHTP-6 S/G	12" X 18"	2000
	LHTP-7 S/G	18" X 24"	3000

**OPTIONAL ACCESSORIES :**

A. Dual display Microprocessor based PID temperature controller & PT – 100 sensor.