



- Useful for Shaking in Research Institutions, Chemical & Agriculture Laboratories.
- The housing made of CRCA Steel - powder coat finished.
- Universal Platform made from S.S. 304 grade to accommodate interchangeable lotus clamps of assorted sizes for different capacity of flasks.
- Shaking Speed can be controlled by Digital RPM Indicator cum Controller between 20 RPM to 250 RPM.
- Permanent Magnet D. C. (PMDC) geared motor for shaking with eccentric orbital drive.
- Shaking Stroke: 25 MM.
- Operates on 230 Volts AC Single phase 50 Hz.

MODEL	UNIVERSAL PLATFORM SIZE	MAXIMUM NO OF CLAMPS			
		100 ML	250 ML	500 ML	1000 ML
LROS-1	350 X 350 MM	16	9	9	4
LROS-2	450 X 450 MM	25	16	16	9
LROS-3	600 X 600 MM	64	36	25	16
LROS-4	750 X 750 MM	81	49	36	25
LROS-5	900 X 900 MM	100	64	49	36

OPTIONAL ACCESSORIES :

A. Digital Timer for Timed shaking process of 999 Minutes.

N. B.: A. S. S. Lotus clamps for holding the flasks in the shakers to be ordered separately.

B. Reciprocal shakers are also available.