



- The complete Soxhlet Extraction Apparatus consists of solvent reservoir extraction assembly with 3, 4 & 6 extraction heaters of various diameters to suit 250 ML flasks to 1 LTR flasks.
 - Outer casing is made of CRCA Steel - finished with powder coated paint.
 - Gap between casing & heater is filled with high grade glass wool insulation.
 - Extraction Heaters made of Fiber glass hand knitted netting having Nichrome heating elements.
 - Heating Elements having individual Regulators & Pilot Lamps to indicate working of the heaters.
 - Control box is away but junctioned with the heating console.
 - Complete with Two retort rods, Cross bar.
 - Operates on 230 Volts AC Single phase 50 Hz.
- N. B.: Glass parts and accessories to be ordered separately.

MODEL	SSEU-3A	SSEU-3B	SSEU-3C	SSEU-4A	SSEU-4B	SSEU-4C	SSEU-6A	SSEU-6B	SSEU-6C
Flask Capacity	250 ML	500 ML	1000ML	250 ML	500 ML	1000ML	250 ML	500 ML	1000ML
Flask Recesses	3	3	3	4	4	4	6	6	6
No of Regulators & Switch	3	3	3	4	4	4	6	6	6
Total Rating	450 W.	600 W.	900 W.	600 W.	800 W.	1200W.	900W.	1200 W.	1800W.