



- Most Advanced Bench top ultrasonic cleaner for surgeons, hospitals & laboratories put an end to the potentially dangerous, inefficient, ineffective hand-scrubbing of delicate expensive instruments by adopting the cleaning technology used world-wide.
- SMPS power source ensures amazing > 90% efficiency, compared to 60% of conventional linear source.
- Intense cleaning action, self tuning, simple to operate.
- S.S. Lid for low noise.
- Presettable 15 minutes digital countdown timer with LED displays.
- Transducer frequency : 33 KHz.
- Supplied with water drain Valve.
- Operates on 230 Volts AC Single phase 50 Hz.

MODEL	TANK DIMENSION (LXBXH) MM	CAPACITY (LITERS)	ULTRASONIC POWER IN WATTS
LUSC - 1	240 x 140 x 65	1.5	50
LUSC - 2	235 x 135 x 100	2.25	70
LUSC - 3	300 x 150 x 100	3.5	100
LUSC - 4	250 x 150 x 150	5.5	150
LUSC - 5	300 x 150 x 150	6.5	200
LUSC - 6	300 x 180 x 150	8.0	200
LUSC - 7	300 x 200 x 150	9.0	250
LUSC - 8	300 x 230 x 150	10.0	250
LUSC - 9	300 x 300 x 150	12.0	300
LUSC - 10	350 x 300 x 200	21.0	400
LUSC - 11	400 x 400 x 150	24.0	500
LUSC - 12	400 x 300 x 250	30.0	500