



- Manesty Type distillation unit designed to provide Pyrogen free distilled water easily.
- Complete unit made of non corrosive Stainless Steel – Argon Welded.
- The lid perfectly rests and gets inside the collar on the head of the tank which avoids loss of steam and jumping of lid.
- Condenser can be opened easily for periodic cleaning.
- The Unit is fitted with immersion type heating element.
- It is provided with iron bracket to be fixed on the wall.

MODEL	CAPACITY LITRS / HR	RATING
LWDU-1	2 LITERS	1500 W.
LWDU-2	4 LITERS	3000 W.
LWDU-3	6 LITERS	4500 W.
LWDU-4	10 LITERS	8000 W.